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
### **DISPOSITIVO DE TROMBECTOMIA pRESET LITE**

*Dispositivo para Neurocirurgia*

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REVISÃO: 02

## English Directions for Use

Legal Manufacturer –  
Phenox Ltd,  
Kamrick Court,  
Ballybrit Business Park,  
Galway H91 XY38  
Ireland.  
Tel: +353 91 740100  
www.phenox.ie

### Content

1 X pRESET LITE Thrombectomy device, consisting of the Stent connected to the pusher wire, covered by the Shrink Tubing of PTFE. The product also has an introducer sheath and a protective cap for the wire.

### Device Description

The pRESET LITE Thrombectomy Device is designed to restore blood flow in the neurovasculature (internal carotid artery (ICA), M1 or M2 segments of the middle cerebral artery (MCA)) by mechanical removal of thrombus in patients experiencing ischemic stroke due to acute vessel occlusion with thrombus.

As per Figure 1:

The pRESET LITE Thrombectomy Device consist of a self-expanding Nitinol structure (1) that carries one radiopaque Platinum/Iridium (Pt/Ir) coil-marker on its proximal end (2) and two markers on its distal end (3). The structure is firmly attached to a pusher wire on its proximal end (4). This pusher wire is covered with a PTFE tubing in order to reduce friction during use. The retrieval device is stored in compressed form in an introducer sheath (not displayed in Figures below).

The pRESET LITE Thrombectomy Device is introduced into the target vessel through a suitable microcatheter (ID 0.0165"") and deployed inside the thrombus or distal to the thrombus. After complete deployment, the slow withdrawal of the retrieval device occurs under aspiration via the guiding catheter.


Figure 1: Representative illustration of the pRESET LITE device

### Indications for Use

The device is indicated for mechanical clot retrieval from intracranial arteries as acute ischemic stroke treatment within 6 hours of symptom onset. Candidates for treatment are patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA), patients who fail IV t-PA therapy and patients who require a supplement treatment to an initiated IV t-PA therapy.

### Contraindications

The device is contraindicated in

- Patients with pre-existing haemorrhage in the brain tissue supplied by the target vessel;
- Patients with known sensitivity to Nickel-Titanium;
- Patients with established infarction in the brain tissue supplied by the target vessel;
- Vessels with a diameter smaller than the lower limit defined in this DFU;
- In tortuous intracranial vessels precluding safe retraction of the device.

## Precautions

- Patients with a Nickel allergy may have a hypersensitivity response to the device. Ensure to weigh the benefit of treatment v's the severity of the "nickel allergy" risks prior to treatment when using pRESET LITE in patients.
- On the strength of past clinical experience, the inherent risks of treatment with pRESET LITE are less in case of an underlying stenosis than in case of a dissection or vasospasm.
- The pRESET LITE can be used up to three (3) times in the same patient. The contact with the patient of all applications must in sum be less than sixty(60) minutes.
- Between the applications, the device must be cleansed of blood and thrombus and inspected for damages. If the device or the introducer sheath is damaged or the device cannot be pulled into the introducer sheath it has to be discarded.
- The device must not be re-sterilized or recycled for other patients because it cannot be reliably cleaned.
- Use before date of expiry.
- Do not use damaged packages.
- Do not deploy the device in an artery smaller than the corresponding minimum vessel diameter according to the table given in the compatibility charts.
- The pRESET LITE is a delicate device and demands cautious treatment. Do not advance or withdraw the device against resistance. Do not rotate the insertion wire. If necessary, remove pRESET LITE via microcatheter.
- The device must only be used by specialized and trained physicians.
- The deployed pRESET LITE must not be retracted through implanted vascular implants.
- Before initial use never push the pRESET LITE out of the introducer sheath (transportation tube) in order to inspect it; pulling the device back would put unnecessary strain on the device.
- Copious flushing of the introducer sheath (transportation tube) is required to remove potentially trapped air.
- Pushing the device beyond the tip of the microcatheter may dissect or perforate the target vessel.
- Pulling or pushing on the insertion wire with inadequate force and torsion of the insertion wire may disconnect the pRESET LITE from the insertion wire. In such cases, the retrieval by a grasp device (Alligator Device) or snare device is recommended.
- In case of excessive retraction resistance, the reintroduction of the microcatheter for splinting the pRESET LITE insertion wire is recommended. If the microcatheter cannot be pushed completely over the device remove both, device and microcatheter as a unit.
- Open the RHV of the guiding catheter completely before pulling the device out.
- Thrombectomy procedures using pRESET LITE is to be carried out in a Clinical operating theatre with fully trained physicians.
- Fluoroscopic visualisation is recommended for confirming correct placement of pRESET LITE during thrombectomy procedures as per step 4 of procedure.
- If vasospasm is suspected in the affected vessel region, all necessary counteractive measures (e.g. medication treatment) should be taken before the pRESET LITE treatment.
- The pRESET LITE device has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment as this is not an implantable device.

## Complications

The following and further complications may occur during or despite the appliance of the pRESET LITE:

- Air-, Catheter-, Plaque-, Thrombus- Emboli;
- Allergic / Toxic Reaction;
- Adverse Response by patient;
- Cerebral Ischemia;
- Extracranial Haemorrhage
- Extracranial / Intracranial Complications during procedure;
- Infection;
- Inflammation;
- Ischemic Stroke/Haemorrhagic Stroke;
- Neurovascular/Neurological complications;
- Severe Harm (permanent disability/death);
- Vasospasm;
- Vessel Perforation / Occlusion / Dissection / Thrombosis / Damage.

### Compatibility

The following pRESET LITE devices and catheters are compatible and are recommended for use in minimum vessel diameters:

Model Number	Device Size	Insertion Wire Length	Microcatheter	Min. vessel Ø	Relaxed Ø x Usable Length
PRE-LT-3-20-B	PRE-LT-3-20	2.0 m	ID 0.0165": - ECHELON 10 / 14 - EXCELSIOR SL 10 - HEADWAY 17	2 mm	3 x 20 mm
PRE-LT-4-20-B	PRE-LT-4-20	2.0 m	ID 0.0165": - ECHELON 10 / 14 - EXCELSIOR SL 10 - HEADWAY 17	2 mm	4 x 20 mm

### Other Accessories for performing the procedure (not supplied):

Microcatheter – refer to compatibility table,  
 Guidewire,  
 Aspiration device (syringe),  
 Saline solution/heparin saline continuous flush set,  
 Rotating Haemostasis valve,  
 Infusion Stand,  
 Femoral Arterial Lock,  
 Guide Catheter,  
 Sheath/long sheath

### Recommended procedure:

1. Preparation for use - remove pRESET LITE thrombectomy device from its package and inspect to ensure no damage to the device.
2. It is important that the device is introduced in sufficiently large target vessels only. Heed the minimum vessel diameter and the catheter compatibility.
3. Introduce a suitable guidewire and microcatheter (see table above) into the target vessel through a Guide Catheter, with careful passage beyond the thrombus. The application of the so-called "road map" technique is recommended. Never catheterize against any resistance! Withdraw the micro guidewire from microcatheter. The careful injection of a small amount (0.2 ml) of contrast medium may be used to verify the proper position of the tip of the microcatheter beyond the thrombus and in the vessel lumen.
4. With a rotating haemostatic valve (RHV) and slow but constant saline flushing connected to the microcatheter, the pRESET LITE will be introduced into the microcatheter. During this process, the RHV will be opened, the introducer sheath with the device inside will pass the valve, the valve will gently be closed and the tube will be flushed back until water exits from the proximal end of the tube. After the tube is flushed with the device inside, the tube will be advanced until it is in close contact with the inner funnel end of the hub of microcatheter. The tube will be firmly kept in this position; the pRESET LITE is advanced from the tube into the microcatheter by gently pushing the insertion wire. The process of device advancement will be continued until about sixty (60) cm of the insertion wire is pushed forward. The tube will be removed. The pRESET LITE is pushed further forward until the marker ("MARKER": distance to distal device tip approximately 125 cm in compressed condition) of the insertion wire reaches the entry of the haemostatic valve. Placing of the device (up to this point) does not require fluoroscopy because the marker ("MARKER") identifies the position to which the system can be advanced without the device tip leaving the microcatheter. However, Fluoroscopy is required to confirm correct positioning of the device (at the target lesion).

**WARNING:** Should major resistance be encountered, remove the device and the microcatheter, start the procedure from the beginning.

5. Deploy the device completely by gentle and slow withdrawal of the microcatheter until the

microcatheter tip is located sufficiently proximal (Figure 2). This is achieved if the distal microcatheter marker, is placed proximal to the proximal x-ray visible pRESET LITE marker. During this process, use continuous fluoroscopy to keep the pRESET LITE in its original position. The instrument is ideally deployed in the thrombus. If the pRESET LITE is longer than the thrombus, the instrument should protrude to the distal end of the thrombus.

**WARNING:** Never push the pRESET LITE beyond the distal tip of the microcatheter, to do so may result in dissection or perforation of the concerning artery.


Figure 2: Device deployment

6. Wait five (5) minutes before starting the withdrawal. The saline flush will be removed from the guiding catheter RHV and will be replaced by a fifty (50) ml syringe with Lure Lock. Two (2) additional syringes of the same size will be prepared and available. Take care to ensure valves are firmly closed to avoid aspiration of air.
7. If applicable, inflate the balloon of the guide catheter according to its DFU.
8. Apply aspiration on the guiding catheter, the unit (microcatheter + pRESET LITE will SLOWLY be pulled back (figure 3) into the guiding catheter. Some resistance may be observed at this point; the more distal the device has been introduced, the more resistance that may be encountered. Always pull gently, never abruptly or in steps!


Figure 3: Revascularization unit Recovery

9. Monitor the device under fluoroscopy as it is travelling through the vessel, under aspiration. As soon as the pRESET LITE has entered the guiding catheter, withdrawal can be accelerated. Pull the device in front of the RHV. Disconnect the RHV and the pRESET LITE from the guiding catheter. The aspiration (with a syringe) is started slowly with the pulling back of the device, this is increased from start to finish and will yield 100 – 150ml of blood. Aspirate for safety reasons another fifty (50) ml of blood from the guiding catheter. As needed close the guiding catheter e.g. using a syringe. If applicable, deflate the balloon of the guide catheter according its DFU
10. Gently inject ten (10) ml of contrast medium into the guiding catheter to see if the above described procedure restored vessel patency.
  - a. If the target vessel is patent, decide if further vessels are occluded and need recanalization. If this is not the case, withdraw the guiding catheter, wait for ten to fifteen (10-15) minutes and re-inject the concerning supra-aortic vessel to confirm permanent patency.
  - b. If the target vessel is not patent, decide upon further action, which may include a repeated pRESET LITE procedure or use of another mechanical device.

11. In case of a repeated or an additional procedure an intact device can be reused in the same patient.

**WARNING:**

1. Do not use the pRESET LITE in more than one patient. Structural integrity and/or function may be impaired by reuse in another patient.
2. Reuse of device in the same patient can be used up to three times. So long as necessary cleaning as detailed in step 12 to 15 must be carried out.

12. Cleaning: Gently rinse the pRESET LITE using sterile water or saline solution; blood or thrombus must not adhere to the device.

13. Inspection: Check the device to exclude any damage.

14. Preparation for new procedure: Push the proximal end of the insertion wire into the corresponding tube (introducer sheath). Then move the tube up to the device. Then pull the device SLOWLY and completely back in the introducer sheath. If the device cannot be pulled into the sheath it has to be discarded.

15. The reapplication follows steps 1 to 11 of this DFU.




**Storage**

Keep away from heat. Store in a cool, dry place.

**Disposal**

The device, retrieved clots, blood and if applicable its packaging must be disposed as biohazard waste.

**Symbols Glossary**

	Prescription only
	This product is sterile and was sterilized using ethylene oxide. Do not use if the <u>sterile packaging is damaged</u> .
	Please read and follow the directions for use before using the product.
	Please read directions for use for cautionary information such as warnings and precautions.
	This product is intended for single use only. Do not reuse!
	Keep product dry and protected from moisture.
	Keep product away from sunlight.
	Do not use the product if package is damaged or opened.
	Do not resterilize the product.
	Product is non-pyrogenic
	Catalogue number
	Batch number (Lot number)
	Expiration date
	Contents
	Legal Manufacturer

**Limits of Liability**

phenox Limited shall not be liable for damage caused by use other than that intended or reuse of the product.

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Phenox Limited